MöllerMedical

**CE 0482**  
Möller Medical GmbH  
Wasserkuppenstraße 29-31 · 36043 Fulda, Germany  
Tel. +49 (0) 661-9 41 95-0 · Fax +49 (0) 661-9 41 95-850  
info@moeller-medical.com · www.moeller-medical.com

## Breast Localisation Needle

Gebrauchsanweisung beachten | Consult instructions for use | Consulter la notice d'utilisation  
Consultare le istruzioni per l'uso | Consultar as instruções de utilização  
Consúltese las instrucciones de uso | Viz návod k použití | Postupujte podľa návodu na použitie

Chargenbezeichnung | Batch code | Code du lot | Codice di lotto  
Código do lote | Código de lote | Šarže | Sériový kód

Bestellnummer | Catalogue number | Référence du catalogue | Riferimento di Catalogo  
Número de catálogo | Número de catálogo | Obj. číslo | Číslo objednávky

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**Indications:** The breast localisation needles serve for the marking and pre-operative identification of breast tissue lesions.

#### Contra-indications:

- This breast localisation needle is exclusively for use by authorised, specialised doctors who are familiar with the potential side effects, typical diagnostic findings, restrictions, indications and contra-indications of markings in the breast tissue.
- The use of the product, in particular on patients with blood coagulation disorders or patients treated with anticoagulants, is at the discretion of the attending doctor.

#### Preparatory measures:

1. Prior to use, make sure that the package is intact and check the use-by date. Only use products available in an undamaged package. The package should be opened under sterile conditions.  
**NOTICE:** Products in an open, damaged or moist package must not be used.
2. Take the breast localisation needle out of the package and then remove the protective cover.
3. Check the the calibrating wire and the cutting edge of the breast localisation needle for damage.  
**NOTICE:** The breast localisation needle with calibrating wire must not be used in case of damage or any other defect which can adversely affect its intended use.
4. **CAUTION:** In order to minimise the risk of deformation of the breast localisation needle, insert the needle into the breast tissue with the utmost care.

**Special features:** All BLM breast localisation needles provide centimetre graduations for monitoring the puncture depth. They differ from each other through the following features:

**BLM-1110:** Stable and also mobile, J-curved, repositionable NiTiNol wire. The length of the J-curved wire is pre-set by means of a chuck. A screw fitting on the proximal end provides for optimal handling and maximum safety when pushing the calibrating wire into the tissue. **NOTICE:** In case of a nickel allergy, the use of the BLM-1110 is at the discretion of the attending doctor.

**BLM-1120:** Smooth, rigid, non-repositionable calibrating wire for optimal tissue penetration. A clip fixed on the skin prevents the dislocation of the calibrating wire. The calibration markings indicate the position of the wire in the breast localisation needle and the activation of the hook. Contrast medium can be injected and liquid can be conveniently aspirated via the open distal end of the breast localisation needle. **NOTICE:** After removing the calibrating wire, make sure that the hook is available.

**BLM-1125:** Same as BLM-1120, but with echogenic tip.

**BLM-1130:** The structured, flexible, non-repositionable calibrating wire provides high flexibility when penetrating the tissue. A clip fixed on the skin prevents the dislocation of the calibrating wire. The calibration markings indicate the position of the wire in the breast localisation needle and the activation of the hook. Contrast medium can be injected and liquid can be conveniently aspirated via the open distal end of the breast localisation needle. **CAUTION:** This calibrating wire must not be cut through; otherwise, there is a risk that the wire unwinds, which can lead to complications and injury. **NOTICE:** After removing the calibrating wire, make sure that the hook is available.

**BLM-1135:** Same as BLM-1130, but with echogenic tip.

**BLM-1140:** Stable but flexible, repositionable NiTiNol wire with a double hook. The length of the wire is preset using a collet. A proximally attached screw fitting ensures optimal handling and maximum safety when advancing the marker wire into the tissue. **CAUTION:** In the event of a nickel allergy, the use of the BLM-1140 is subject to the attending doctor's discretion.

#### General remarks:

- The product must be kept away from light and stored in a dry and cool place.
- The attending doctor is responsible for the proper use and execution of the examination using this product (these operating instructions only serve as a general guideline). A calibrating wire may only be removed by an experienced surgeon.
- This breast localisation needle is a sterile single-use product intended for use on one patient. Resterilisation and reuse are not permitted. The reuse of single-use products involves a potential risk of infection for the patient or operator. Contaminated products can cause damage to health, disease or even death of the patient. Cleaning, disinfection or sterilisation can adversely affect essential material properties and product parameters and thus result in product failure.
- The product must be disposed of in accordance with hospital regulations and applicable laws.
- Möller Medical GmbH expressly disclaims any liability for damage caused through improper use of the products or use by unauthorised personnel.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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