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Description The Bone Access Needle is a sterile packaged kit for gaining access to bone.

Intended use The bone access needle is a sterile-packed kit used in procedures where access to bone is required. The bone access needle can also be used for providing access for appropriate guide wires, for obtaining bone marrow using the aspiration technique and for inserting therapeutic materials, including bone cements. Refer to the corresponding instruction for use regarding indications, contraindications, compatibility, use, precautions, warnings and side effects when using the bone access needle in connection with other instruments and/or therapeutic materials, including bone cements. For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended.

Indications The indications for the Bone Access Needle are defined by the procedure for which access to bone is required and can be found in the corresponding instruction for use.

The use of the Bone Access Needle may be indicated by the necessity to extract bone marrow cells, for instance, for histological examination, which is at the discretion of the attending physician, or by the procedure for which access to bone is required and can be found in the corresponding instruction for use.

Contraindications This device is not to be used in patients with coagulation disorders or infections.

Patient Target Group The device is applicable to all patients who, at the discretion of the physician, may be treated in compliance with the intended use, indications, contraindications, warnings and precautions.

Warnings and precautions Since the device is used in technically complex procedures, the physicians should be familiar with the handling and use of the device, the other instrumentation and the procedure.

Before carrying out the procedure, the physician must ensure that the size of the available device is appropriate for the patient.

It is also necessary to check whether the patient has allergies to materials (e.g. nickel) contained in the device. Should this be the case, the use of the device is at the discretion of the attending physician.

The safe use of the device cannot be guaranteed if used in conjunction with MRI, which is associated with certain risks, including heating or migration of the device and artefacts on the MRI image.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (only for European Union). Other countries should follow local regulations or procedures.

Side effects No side effects are known to date.

Materials The device contains the following materials:

- Stainless steel: Needles and mandrins
- Polybutylene terephthalate and Polyamide 12: Handles and other plastic components

Device variants and content The needles are available with blue (8G), violet (11G) or orange (13G) handles. Each device contains a needle (with Luer-lock) and a mandrin.

Mode of Use The correct placement of all instruments and needles must be performed with the aid of imaging procedures.

The Bone Access Needle can be inserted into the bone with gentle taps of the hammer.

CAUTION: Strong hammer blows should be avoided, even if sclerotic or thicker cortical bone is involved. Under certain circumstances, there might be a risk of splitting the bone.

Sterility The device is supplied sterile. Before opening, check that the packaging is intact. Remove items from the packaging while maintaining aseptic conditions. Do not use if the packaging seal is broken or if the sterile barrier has been breached.

Single-use The device is intended for single-use only and may not under any circumstances be reused. The reuse of a single-use device poses an infection risk to patients and users. Contamination of the device can lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization can adversely affect important material and design properties and device function.

Storage The device should be stored in a clean, dry environment and protected from direct sunlight. Do not use after the expiration date on the label.

Disposal Disposal takes place in accordance with the hospital ordinance or the applicable laws. Due to the risk of infection and injury, care must be taken to avoid all contact with the tips or sharp edges of the instruments. This is especially important for contaminated instruments.

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