Möller Medical GmbH
Wasserkuppenstraße 29-31
36043 Fulda, Germany
Tel.: +49 (0) 661-9 41 95-0
Fax: +49 (0) 661-9 41 95-850
info@moeller-medical.com
www.moeller-medical.com

de

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en

Instructions for use

fr

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it

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es

Folleto de instrucciones

pt

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nl

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da

Brugsanvisning

sv

Bruksanvisning

no

Bruksanvisning

fi

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tr

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el

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Möller Medical Blue

en Instructions for use for reusable, automatic biopsy guns

Möller Medical Blue
 (REF.: RBG-1000-10-1000)

Explanation of safety symbols used

Important information is accompanied by visual symbols in these instructions for use. This information is essential for avoiding risks to the patient and operating personnel and for avoiding biopsy gun damage and malfunctions.

Purpose

The Möller Medical Blue needles and systems are used for taking and suctioning tissue samples. The various wall thicknesses, tips, dimensions and cutting techniques of the different core needle types and biopsy systems all influence the quality of the biopsy specimen. Product applications are soft tissue samples, e.g., kidneys, liver, prostate, lymph nodes, chest and bones/bone marrow.

Indications:

Möller Medical Blue is a reusable, automatic biopsy gun that is intended exclusively for use with DNG-1020 / full-sCore® DNG-2020 disposable needles for core needle biopsies of soft tissue.

Contraindications:

- The use of these biopsy guns is restricted to registered, specialist doctors familiar with the possible side effects, typical findings, restrictions, indications and contraindications of core needle biopsies.
- Performance of a core needle biopsy using these biopsy guns, especially in patients with a nickel allergy, blood-clotting disorders or patients currently being treated with anticoagulants, is subject to the attending doctor's discretion.

Preparation of the Blue biopsy gun (see Fig. 1):

- Remove the biopsy gun from its packaging.
-  The biopsy gun is not sterile upon delivery and must be cleaned and sterilised before each use in accordance with the reprocessing details given in these instructions.
- Test the biopsy gun before use: Cock the biopsy gun without the needle inserted by pulling the lever (J) backwards twice. Check the function of the biopsy gun by first releasing the safety catch (C) and then pressing the front or rear trigger (D).


Figure 1: Möller Medical Blue

-   Never check the biopsy gun with the needle inserted. This could result in damage to the needle and/or injury to the patient or user. Use the biopsy gun only with the cover closed.

Use:

1. Before performing a core needle biopsy, check that the packaging is undamaged and that the expiry date of the DNG-1020 / full-sCore® DNG-2020 disposable biopsy needle has not elapsed. Use only products with undamaged packaging. Apply aseptic measures when opening the packaging.

 Do not use products if their packaging is already open, damaged or damp.

2. When using the DNG-1020 or full-sCore® DNG-2020, note the relevant instructions for the disposable biopsy needles for use with biopsy guns.

 If you identify any damage or other defects that affect the correct functioning of the disposable biopsy needles, they must not be used.

3. Inserting a biopsy needle

- 3.1 **Insert the DNG-1020 disposable biopsy needles (see Fig. 2 and 3) as follows:** Hold the gun in your right hand and open the cover (A) of the automatic biopsy gun from left to right. Insert the DNG-1020 disposable biopsy needle. Make sure that the plastic supports of the DNG-1020 (E) are securely attached on the locking pins (F, G). Now pull out the spacer (B) as you reclose the cover (A) – the cover must click into place. Remove the needle's protective sheath.


Figure 2: Inserting the DNG-1020


Figure 3: Removing the DNG-1020 spacer

3.2 Insert the full-sCore® DNG-2020 (see Fig. 4 and 5) as follows:

Hold the biopsy gun in your right hand and open the cover (A) of the automatic biopsy gun from left to right. Insert the full-sCore® DNG-2020. Please ensure that the plastic support (E) of the full-sCore® DNG-2020 is positioned carefully on the front locking pin (F) while the rear part of the plastic support (H) of the full-sCore® DNG-2020 fits flush against the centre block (I); the rear locking pin remains free. Now pull out the spacer (B) (see Fig. 5) as you reclose the cover (A) – the cover must click into place. Remove the needle's protective sheath.

 To avoid infections, always work in aseptic conditions, including disinfecting the skin, wearing sterile gloves and using sterile drapes as necessary.


Figure 4: Inserting the DNG-2020


Figure 5: Removing the DNG-2020 spacer

4.  The biopsy gun is made ready to fire by applying medium force to pull back the lever (J) twice. To do this, hold the biopsy gun in your hand. Grasp the bottom part of the lever (J) with the index finger of your other hand, and place your thumb on the top part of the lever (J). The device is now automatically secured. The safety catch is in the LOCKED position.
5. The lock can be released by pushing the safety catch (C) to the UNLOCKED position.
6. Insert the tip of the needle under imaging guidance into the patient up to just before the lesion.
7. Press the front (D) or rear trigger (D). This activates a spring mechanism and the needle for taking the biopsy is released.

8. Remove the needle from the patient. To release the sample from the core needle, the DNG-1020 disposable or full-sCore® DNG-2020 biopsy needle need not be removed from the biopsy gun. You can release and remove the tissue sample by pressing the lever (J) once. In the same step, the safety catch (C) is automatically moved to the LOCKED position.
9. Remove the contaminated disposable biopsy needle from the biopsy gun by pressing the lever (J) again to activate the biopsy gun. Once the trigger (D) is pressed, the needle is in its basic position. You can now remove it and dispose of it (see Fig. 6).


Figure 6

 Risk of infection from contaminated needles!

Troubleshooting

Problem	Possible cause	Remedy
Biopsy gun is difficult to cock or cannot be fully cocked	<ul style="list-style-type: none"> - Improper handling - Lever has not been pulled back completely - Safety catch is in the wrong position - Biopsy gun is soiled - Device is faulty 	<ul style="list-style-type: none"> - Correct handling - Always pull back the lever fully - Check the safety catch - Clean and lubricate the biopsy gun - Return to manufacturer
Biopsy gun does not release or is difficult to release	<ul style="list-style-type: none"> - Biopsy gun is not fully cocked - Biopsy gun is not unlocked - Biopsy gun is soiled - Lubrication is insufficient - Device is faulty 	<ul style="list-style-type: none"> - Check whether the biopsy gun is fully cocked - Check whether the biopsy gun is unlocked - Clean and lubricate the biopsy gun - Return to manufacturer
Biopsy gun does not release fully	<ul style="list-style-type: none"> - Biopsy gun is soiled - Biopsy gun is not sufficiently lubricated - Device is faulty 	<ul style="list-style-type: none"> - Clean and lubricate the biopsy gun - Return to manufacturer

Möller Medical recommends cleaning and lubricating the Möller Medical Blue after each use.

All reprocessing instruction and information as given in the IFU.

General information:

- Use the biopsy gun only with the cover closed.
- The biopsy gun can be triggered by severe impacts when it is cocked (e.g. falling on the floor). When setting the biopsy gun down, make sure the surface is flat, secure and stable.
- The attending doctor is responsible for the proper use and performance of the examination with these products (these instructions for use are only general instructions).
- Möller Medical GmbH accepts no liability for damage due to improper use of the products, the use of disposable biopsy needles from third-party manufacturers or the use of the products by unauthorised personnel.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Contact the Manufacturer: Möller Medical GmbH
Wasserkuppenstraße 29-31
36043 Fulda, Germany
Tel.: +49 (0) 661-9 41 95-0
Fax: +49 (0) 661-9 41 95-850
info@moeller-medical.com
www.moeller-medical.com

Temperature range for use: 0 °C to 50 °C

 The DNG-1020 and full-sCore® DNG-2020 disposable biopsy needles are designed for **single-use only** and must not be reprocessed!

 **For USA only:** When used for breast biopsy, the product is for diagnosis only. The extent of histological abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histological abnormality; e.g., malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal of using standard surgical procedures.

 Attention: Under US federal law, this device may be only sold to a physician or ordered by a physician.

Möller Medical GmbH
Wasserkuppenstr. 29-31 · 36043 Fulda, Allemagne
Tel. +49 (0) 661-9 41 95-0 · Fax +49 (0) 661-9 41 95-850
info@moeller-medical.com · www.moeller-medical.com

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