MöllerMedical

CE 0482

Möller Medical GmbH

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Instructions for Use icotec Instrument Cement Cannula

Device in Scope

REF.	Device Name
42-420	icotec Instrument Cement Cannula

Description The icotec Instrument Cement Cannula is a sterile packed cannula designed for the augmentation of fenestrated pedicle screws.

Purpose The icotec Instrument Cement Cannula is intended for cement augmentation in combination with the fenestrated pedicle screws of the icotec Pedicle Systems. For further information, please refer to the respective Surgical Technique manual and instructions for use of the icotec Pedicle Systems.

Indications The indications for the icotec Instrument Cement Cannula are defined by the icotec Pedicle Systems and can be found in the respective instructions for use.

Contraindications This device is not to be used in patients with coagulation disorders or infections. Further contraindications are defined by the icotec Pedicle Systems and can be found in the respective instructions for use.

Patient Target Group The patient target group is defined by the icotec Pedicle Systems and can be found in the respective instructions for use.

Warnings and precautions Since the device is used in technically complex procedures, the physicians should be familiar with the handling and use of the device, the other instrumentation and the procedure.

It is also necessary to check whether the patient has allergies to materials (e.g. nickel) contained in the device. Should this be the case, the use of the device is at the discretion of the attending physician.

The safe use of the device cannot be guaranteed if used in conjunction with MRI, which is associated with certain risks, including heating or migration of the device and artefacts on the MRI image.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (only for European Union). Other countries should follow local regulations or procedures.

Side effects The side effects are defined by the icotec Pedicle Systems and can be found in the respective instructions for use.

Materials The device contains the following materials:

- Stainless steel: Cannula and Plunger
- Polybutylene terephthalate: Handles and other plastic components

Device content The device contains the following items:

- Cannula with shaped tip (1 unit)
- Plunger (1 unit)

Mode of use For further information on the use of the icotec Instrument Cement Cannula, please refer to the respective Surgical Technique manual and instructions for use for the icotec Pedicle Systems.

Sterility The device is supplied sterile. Before opening, check that the packaging is intact. Remove items from the packaging while maintaining aseptic conditions. Do not use if the packaging seal is broken or if the sterile barrier has been breached.

Single-Use Device The device is intended for single use only and may not under any circumstances be reused. The reuse of single-use devices poses an infection risk to patients and users. Contamination of the device can lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization can adversely affect important material and design properties and device function.

Storage The device should be stored in a clean, dry environment and protected from direct sunlight. Do not use after the expiration date on the label.

Disposal Disposal takes place in accordance with the hospital ordinance or the applicable laws. Due to the risk of infection and injury, care must be taken to avoid all contact with the tips or sharp edges of the instruments. This is especially important for contaminated instruments.

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