# CE 0482

Möller Medical GmbH

Wasserkuppenstraße 29-31 · 36043 Fulda, Germany  
Tel.: +49 (0) 661-9 41 95-0 · Fax: +49 (0) 661-9 41 95-850  
info@moeller-medical.com · www.moeller-medical.com

Distributed by

**DePuy Synthes**

COMPANIES OF Johnson & Johnson

- de** Gebrauchsanweisung
- en** Instruction leaflet
- fr** Notice d'instructions
- it** Foglietto di istruzioni
- es** Folleto de instrucciones
- pt** Folheto de instruções
- hu** Használati utasítás
- cs** Návod k použití
- sk** Návod na použitie

[www.moeller-medical.com/ifu-ga](http://www.moeller-medical.com/ifu-ga)

- bg** Инструкция за употреба
- da** Brugsvejledning
- el** Οδηγίες χρήσης
- et** Kasutusjuhend
- fi** Käyttöohje
- hr** Upute za uporabu
- lt** Naudojimo instrukcija
- lv** Lietošanas instrukcija
- nl** Gebruiksaanwijzing
- no** Bruksanvisning
- pl** Instrukcja obsługi
- ro** Instrucțiuni de utilizare
- si** Navodila za uporabo
- sv** Bruksanvisning
- tr** Kullanım kılavuzu

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## Instructions for Use TRAUMACEM™ V+ Injection Cannula

### Device in scope

03.702.120S	TRAUMACEM™ V+ Injection Cannula, for PFNA System
03.702.121S	TRAUMACEM™ V+ Injection Cannula, for TFNA System

**Description** TRAUMACEM™ V+ Injection Cannulas for the TFNA and PFNA Systems are sterile packed devices for the augmentation of suitable DePuy Synthes implants.

**Intended use** The TRAUMACEM™ V+ Injection Cannula is intended for the augmentation of suitable DePuy Synthes implants with PMMA<sup>1</sup>-based bone cement into the proximal femur region. The TRAUMACEM™ V+ Injection Cannula may only be used with the TRAUMACEM™ V+ Injectable Bone Cement and the TRAUMACEM™ V+ Syringe Kit.

Refer to the corresponding instruction for use regarding indications, contraindications, compatibility, use, precautions, warnings and side effects of the bone cement used in conjunction with the TRAUMACEM™ V+ Injection Cannula.

For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended.

**Indications** The indications for the augmentation procedure using the TRAUMACEM™ V+ Injection Cannula are defined by the PMMA-based bone cement and can be found in the respective instructions for use.

**Contraindications** This device is not to be used in patients with coagulation disorders or infections.

Further contraindications are defined by the PMMA-based bone cement and can be found in the respective instructions for use.

**Patient Target Group** The patient target group is defined by the PMMA-based bone cement and can be found in the respective instructions for use.

**Warnings and precautions** Since the device is used in technically complex procedures, the physicians should be familiar with the handling and use of the device, the other instrumentation and the procedure.

Before carrying out the procedure, the physician must ensure that the size of the available device is appropriate for the patient.

It is also necessary to check whether the patient has allergies to materials (e.g. nickel) contained in the device. Should this be the case, the use of the device is at the discretion of the attending physician.

<sup>1</sup>PMMA = Polymethyl Methacrylate

The safe use of the device cannot be guaranteed if used in conjunction with MRI, which is associated with certain risks, including heating or migration of the device and artefacts on the MRI image.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (only for European Union). Other countries should follow local regulations or procedures.

**Side effects** The side effects are defined by the PMMA-based bone cement and can be found in the respective instructions for use.

**Materials** The device contains the following materials:

- Stainless steel: Cannula including sleeve and plunger
- Polybutylene terephthalate: Handles and other plastic components

**Device variants and content** Each device contains a side-opening cannula (with Luer-lock) including a sleeve and a plunger.

**Mode of Use** For further information on the use of the TRAUMACEM™ V+ Injection Cannula please refer to the respective DePuy Synthes surgical technique and instruction for use for the Trauma Augmentation systems.

**Sterility** The device is supplied sterile. Before opening, check that the packaging is intact. Remove items from the packaging while maintaining aseptic conditions. Do not use if the packaging seal is broken or if the sterile barrier has been breached.

**Single-use** The device is intended for single-use only and may not under any circumstances be reused. The reuse of a single-use device poses an infection risk to patients and users. Contamination of the device can lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization can adversely affect important material and design properties and device function.

**Storage** The device should be stored in a clean, dry environment and protected from direct sunlight. Do not use after the expiration date on the label.

**Disposal** Disposal takes place in accordance with the hospital ordinance or the applicable laws. Due to the risk of infection and injury, care must be taken to avoid all contact with the tips or sharp edges of the instruments. This is especially important for contaminated instruments.

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