Möller Medical GmbH
Wasserkuppenstraße 29-31 · 36043 Fulda, Germany
Tel. +49 (0) 661-9 41 95-0 · Fax +49 (0) 661-9 41 95-850
info@moeller-medical.com · www.moeller-medical.com

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en Instructions for Use Biopsy Needle

Device in scope

| | | |
|-------------|--------------------------------------|--------------|
| 03.702.250S | Biopsy needle for Bone Access Needle | 8 G, 140 mm |
| 03.702.251S | Biopsy needle for Bone Access Needle | 11 G, 100 mm |
| 03.702.252S | Biopsy needle for Bone Access Needle | 11 G, 150 mm |

Description The Biopsy Needle is a sterile packaged kit for taking bone biopsies.

Intended use The Biopsy Needle is intended for taking bone biopsies and to be used in conjunction with the Möller Medical Bone Access Needle.

For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended.

Indications The use of the Biopsy Needle is indicated in cases where, during a procedure using the Möller Medical Bone Access Needle, a bone biopsy is required for histological examination.

Contraindications This device is not to be used in patients with coagulation disorders or infections.

Patient Target Group The device is applicable to all patients who, at the discretion of the physician, may be treated in compliance with the intended use, indication, contraindications, warnings and precautions.

Warnings and precautions Since the device is used in technically complex procedures, the physicians should be familiar with the handling and use of the device, the other instrumentation and the procedure.

Before carrying out the procedure, the physician must ensure that the size of the available device is appropriate for the patient. If using the device together with other instrumentation, the physician must ensure the correct chosen sizes to guarantee the compatibility.

It is also necessary to check whether the patient has allergies to materials (e.g. nickel) contained in the device. Should this be the case, the use of the device is at the discretion of the attending physician.

The safe use of the device cannot be guaranteed if used in conjunction with MRI, which is associated with certain risks, including heating or migration of the device and artefacts on the MRI image.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (only for European Union). Other countries should follow local regulations or procedures.

Side effects No side effects are known to date.

Materials The device contains the following materials:

- Stainless steel: Needles and styles
- Polybutylene terephthalate: Handles and other plastic components
- High-density Polyethylene: Luer cap

Device variants and content The needles are color coded in blue (for Bone Access Needle 8G) or violet (for Bone Access Needle 10G). Each device contains a biopsy needle (with Luer-lock), a stylet and a Luer cap. The compatibility of the Biopsy Needle with the corresponding Bone Access Needle is shown in the table below.

| Variant | | | Corresponding Product |
|--------------------------------------|--------------|------------------|--------------------------------------|
| Biopsy Needle for Bone Access Needle | 8 G, 140 mm | REF. 03.702.250S | REF. 03.702.240S |
| Biopsy Needle for Bone Access Needle | 11 G, 100 mm | REF. 03.702.251S | REF. 03.702.241S REF. 03.702.243S |
| Biopsy Needle for Bone Access Needle | 11 G, 150 mm | REF. 03.702.252S | REF. 03.702.244S REF. 03.702.245S |

Mode of Use The correct placement of all instruments and needles must be performed with the aid of imaging procedures.

- Insert the Biopsy Needle into the placed needle of the Bone Access Needle and lock it into the handle.
- Slowly advance the needle assembled with the biopsy needle to the final position in the bone. During advancement of the assembly tissue is captured in the biopsy needle.

NOTE: The assembly can be inserted into the bone with gentle taps of the hammer if the Luer cap is screwed onto the biopsy needle.

- Rotate the assembly one full turn.

NOTE: This will help to loosen the bone biopsy.

- Attach a standard syringe (not included) to the Luer-lock of the biopsy needle and slowly draw the syringe to create a vacuum.

- Remove the biopsy needle with the syringe attached from the needle of the Bone Access Needle.

- Remove the syringe from the biopsy needle and use the stylet to push the collected tissue out of the biopsy needle.

Sterility The device is supplied sterile. Before opening, check that the packaging is intact. Remove items from the packaging while maintaining aseptic conditions. Do not use if the packaging seal is broken or if the sterile barrier has been breached.

Single-use The device is intended for single-use only and may not under any circumstances be reused. The reuse of a single-use device poses an infection risk to patients and users. Contamination of the device can lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization can adversely affect important material and design properties and device function.

Storage The device should be stored in a clean, dry environment and protected from direct sunlight. Do not use after the expiration date on the label.

Disposal Disposal takes place in accordance with the hospital ordinance or the applicable laws. Due to the risk of infection and injury, care must be taken to avoid all contact with the tips or sharp edges of the instruments. This is especially important for contaminated instruments.

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