
SBG-2010 Disposable Semi-Automatic Biopsy Gun

i Gebrauchsanweisung beachten | Consult instructions for use | Consulter la notice d'utilisation
 Consultare le istruzioni per l'uso | Consultar as instruções de utilização
 Consultense las instrucciones de uso | Viz návod k použití | Postupujte podľa návodu na použití

LOT Chargenbezeichnung | Batch code | Code du lot | Codice di lotto
 Código do lote | Código de lote | Šarže | Sériový kód

REF Artikelnummer | Catalogue number | Référence du catalogue | Riferimento di Catalogo
 Número de catálogo | Número de catálogo | Obj. číslo | Číslo objednávky

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en Operating instructions for disposable semiautomatic biopsy guns SBG-2010

Indication for Use:

The Möller Medical Biopsy Needles and Systems are used to obtain core samples and aspiration from tissue by biopsy. The various wall thicknesses, diameters, tip geometries and cutting techniques of the different hollow needle types and biopsy systems affect the quality of the biopsy specimen. Application areas of the products are soft tissue such as kidney, liver, prostate, lymph nodes, breast and bone resp. bone marrow.

Indications for the SBG: The disposable semiautomatic biopsy guns are intended for biopsies of soft tissue.

Contra-indications:

- This disposable semiautomatic biopsy gun is exclusively for use by authorised, specialised physicians who are familiar with the potential side effects, typical diagnostic findings, restrictions, indications and contra-indications of hollow needle biopsies.
- The use of the product, in particular on patients with a nickel allergy, blood coagulation disorders or patients treated with anticoagulants, is at the discretion of the attending physician.

Instructions for use:

1. Prior to use, make sure that the package is intact and check the use-by date. Only use products available in an undamaged package. The package should be opened under sterile conditions. **NOTICE:** Products in an open, damaged or moist package must not be used.
2. Take the disposable semiautomatic biopsy gun out of the package and then remove the protective cover.
3. Check the stiletto tip (A) and the needle's cutting edge (B) for damage. **NOTICE:** The disposable semiautomatic biopsy gun must not be used if damaged or defective in any way that could adversely affect its intended use or if the sample kerf of the stiletto protrudes when the disposable semiautomatic biopsy gun is loaded.
4. To prevent infection, make sure to work under sterile conditions, including skin disinfection, sterile covering if applicable and sterile gloves.
5. Place two fingers in the rings and pull the rear knob with the other hand until the latching click of the spring is heard. The needle can be loaded up to the first or the second click. While the first click indicates needle loading with a 10 mm sampling cavity (notch) the second click ensures a 20 mm sampling cavity (notch).
6. After checking that the needle has been loaded in the desired position, insert it in the site selected for the biopsy. Keep the instrument position and the needle orientation unaltered. **NOTICE:** Take care not to push the rear knob until the needle is in the desired position.
7. When the needle is in place, push the rear knob forward until resistance is noticeable.
8. Firmly push the rear knob once more until the final latching click is heard in order to take the material to be analysed.
9. Then gently remove the needle from the patient tissues.
10. To recover the fragment collected for analysis, repeat the loading operation by pulling the rear knob until the latching click of the spring at 10 or 20 mm is heard.
11. Then slowly push the knob for only ¼ of its travel, that way the needle is extracted and the biopsy material can be transferred onto the slide.

General remarks:

- The product must be kept away from light and stored in a dry and cool place.
- The attending physician is responsible for the proper use and execution of the examination using this product (these operating instructions only serve as a general guideline).
- This disposable semiautomatic biopsy gun is a sterile single-use product intended for use on one patient. Resterilisation and reuse are not permitted. The reuse of single-use products involves a potential risk of infection for the patient or operator. Contaminated products can cause damage to health, disease or even death of the patient. Cleaning, disinfection or sterilisation can adversely affect essential material properties and product parameters and thus result in product failure.
- It is also recommended not to use the device if particularly risky manoeuvres are necessary where the needle will be subjected to high torsional or bending forces which might cause it to snap at the height of the grip.
- The product must be disposed of in accordance with hospital regulations and applicable laws.
- Möller Medical GmbH expressly disclaims any liability for damage caused through improper use of the products or use by unauthorised personnel.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

For USA only: When used for breast biopsy, the product is for diagnosis only. The extent of histological abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histological abnormality; e.g., malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal of using standard surgical procedures.

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