CE 0482

Möller Medical GmbH

Wasserkuppenstraße 29-31 · 36043 Fulda, Germany
Tel.: +49 (0) 661-9 41 95-0 · Fax: +49 (0) 661-9 41 95-850
info@moeller-medical.com · www.moeller-medical.com

Distributed by

DePuy Synthes

COMPANIES OF Johnson & Johnson

- de Gebrauchsanweisung
- en Instruction leaflet
- fr Notice d'instructions
- it Foglietto di istruzioni
- es Folleto de instrucciones
- pt Folheto de instruções
- hu Használati utasítás
- cs Návod k použití
- sk Návod na použitie

www.moeller-medical.com/ifu-ga

- bg Инструкция за употреба
- da Brugsvejledning
- el Οδηγίες χρήσης
- et Kasutusjuhend
- fi Käyttöohje
- hr Upute za uporabu
- lt Naudojimo instrukcija
- lv Lietošanas instrukcija
- nl Gebruiksaanwijzing
- no Bruksanvisning
- pl Instrukcja obsługi
- ro Instrucțiuni de utilizare
- si Navodila za uporabo
- sv Bruksanvisning
- tr Kullanım kılavuzu

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Instructions for Use Biopsy Kit

Device in scope

03.702.222S	Biopsy Kit for 8G Vertebroplasty Needle Kit
03.702.223S	Biopsy Kit for 10G Vertebroplasty Needle Kit

Description The Biopsy Kit is a sterile packaged kit for taking bone biopsies from the vertebral body.

Intended use The Biopsy Kit is intended for taking vertebral body bone biopsies during Vertebroplasty procedures and to be used in conjunction with the Möller Medical Vertebroplasty Needle Kit of the Vertecem Vertebroplasty System / Vertecem V+ System¹. For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended.

¹ Vertecem Vertebroplasty System is only available for sale in the United States. Vertecem V+ System is only available for sale outside the United States.

Indications The use of the Biopsy Kit is indicated in cases where, during a Vertebroplasty procedure using the Möller Medical Vertebroplasty Needle Kit, a bone biopsy is required for histological examination.

Contraindications This device is not to be used in patients with coagulation disorders or infections.

Patient Target Group The device is applicable to all patients who, at the discretion of the physician, may be treated in compliance with the intended use, indications, contraindications, warnings and precautions.

Warnings and precautions Since the device is used in technically complex procedures, the physicians should be familiar with the handling and use of the device, the other instrumentation and the procedure.

Before carrying out the procedure, the physician must ensure that the size of the available device is appropriate for the patient. If using the device together with other instrumentation, the physician must ensure the correct chosen sizes to guarantee the compatibility. It is also necessary to check whether the patient has allergies to materials (e.g. nickel) contained in the device. Should this be the case, the use of the device is at the discretion of the attending physician.

The safe use of the device cannot be guaranteed if used in conjunction with MRI, which is associated with certain risks, including heating or migration of the device and artefacts on the MRI image.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (only for European Union). Other countries should follow local regulations or procedures.

Side effects No side effects are known to date.

Materials The device contains the following materials:

- Stainless steel: Needles and stylets
- Polybutylene terephthalate: Handles and other plastic components

Device variants and content The needles are color coded in blue (for use with Vertebroplasty Needle Kit 8G, Diamond Tip) or yellow (for use with Vertebroplasty Needle Kit 10G, Diamond and Beveled Tip). Each device contains a biopsy needle (with Luer-lock) and a stylet. The compatibility of the Biopsy Kits with the corresponding Vertebroplasty Needle Kits is shown in the table below.

Variant		Corresponding Product
Biopsy Kit for 8G Vertebroplasty Needle Kit	REF. 03.702.222S	REF. 03.702.216S
Biopsy Kit for 10G Vertebroplasty Needle Kit	REF. 03.702.223S	REF. 03.702.218S REF. 03.702.219S

Mode of Use The correct placement of all instruments and needles must be performed with the aid of imaging procedures. Anterior/posterior and lateral views must be used.

CAUTION: When placing the instruments and needles it is essential to ensure that the anterior wall of each vertebral body is not penetrated.

Please refer to the Vertecem Vertebroplasty System / Vertecem V+ System Technique Guide for details.

- When the inserted needle assembly of the Vertebroplasty Needle Kit has passed the posterior wall,
 - For Vertebroplasty Needle Kit 8G and 10G, Diamond Tip (REF. 03.702.216S and 03.702.218S): remove the cannulated trocar and the guide wire.
 - For Vertebroplasty Needle Kit 10G, Beveled Tip (REF. 03.702.219S): remove the beveled trocar.
- Insert the biopsy needle into the needle of the Vertebroplasty Needle Kit and lock it into the handle.
- Slowly advance the needle assembled with the biopsy needle to the final position in the anterior half of the vertebral body. During advancement of the assembly tissue is captured in the biopsy needle.

CAUTION: Do not use a hammer and avoid perforating the anterior wall.

- Rotate the assembly one full turn.

NOTE: This will help to loosen the biopsy.

- Attach a standard syringe (not included) to the Luer-lock of the biopsy needle and slowly draw the syringe to create a vacuum.
- Remove the biopsy needle with the attached syringe from the needle of the Vertebroplasty Needle Kit.
- Remove the syringe from the biopsy needle and use the stylet to push the collected tissue out of the biopsy needle.

Sterility The device is supplied sterile. Before opening, check that the packaging is intact. Remove items from the packaging while maintaining aseptic conditions. Do not use if the packaging seal is broken or if the sterile barrier has been breached.

Single-use The device is intended for single-use only and may not under any circumstances be reused. The reuse of a single-use device poses an infection risk to patients and users. Contamination of the device can lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization can adversely affect important material and design properties and device function.

Storage The device should be stored in a clean, dry environment and protected from direct sunlight. Do not use after the expiration date on the label.

Disposal Disposal takes place in accordance with the hospital ordinance or the applicable laws. Due to the risk of infection and injury, care must be taken to avoid all contact with the tips or sharp edges of the instruments. This is especially important for contaminated instruments.

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