
# BOW-1064 Bone Marrow Biopsy System

**EN** **Operating instructions for BOW-1064 bone marrow biopsy systems**

**Intended Use:** The Möller Medical Biopsy Needles and Systems are used to obtain core samples and aspiration from tissue by biopsy. The various wall thicknesses, diameters, tip geometries and cutting techniques of the different hollow needle types and biopsy systems affect the quality of the biopsy specimen.

**Indications:** The BOW-1064 bone marrow biopsy systems are intended for the atraumatic taking of specimens from and inside bones. Application areas of the products are bone resp. bone marrow.

**Contra-indications:** The use of the product, in particular on patients with blood coagulation disorders or patients treated with anticoagulants, is at the discretion of the attending doctor.

**Patient Target Group:** The device is applicable to all patients who, at the discretion of the physician, may be treated in compliance with the intended use, indications, contraindications, warnings and precautions.

**Preparatory measures:**

- This bone marrow biopsy system is exclusively for use by authorised, specialised doctors who are familiar with the potential side effects, typical diagnostic findings, restrictions, indications and contra-indications of bone marrow biopsies.
- Prior to use, make sure that the package is intact and check the use-by date. Only use products available in an undamaged package. The package should be opened under sterile conditions. **NOTICE:** Products in an open, damaged or moist package must not be used.
- Take the products out of the package and then remove the protective covers.
- Check the mandrin tip (5) and the needle's cutting edge (2) for damage and make sure that the mandrin (6) slides freely within the biopsy needle (1). **NOTICE:** Bone marrow biopsy systems must not be used in case of damage or any other defect which can adversely affect their intended use.
- CAUTION:** Do not exert excessive force during application in order to prevent deformation or breakage of the biopsy needle causing inadvertent injury to the patient.

**Intended use, general information:**

- Correct patient positioning.
- Definition of the entry point and marking if required.
- To prevent infection, make sure that you work under sterile conditions, including skin disinfection, sterile covering and sterile gloves.
- Local anaesthesia of the puncture point and the particularly sensitive periosteum.
- If required, perform a stab incision using the scalpel in order to facilitate skin penetration.
- Hold the biopsy needle by the handle (3), place your index finger on the needle and extend the finger towards the distal end. This enables safe handling and control over the maximum penetration depth.
- Using the biopsy needle, exert rotating pressure in order to form a puncture channel which first extends to the periosteum (resistance is noticeable), and then continue with the utmost care until you have overcome the compacta (the resistance suddenly decreases).
- After entering the bone marrow cavity, remove the mandrin from the biopsy needle by rotating the locking mechanism 90°.

Points 1. to 8. apply to both aspiration and biopsy. All further steps are explained in detail in the Intended use ASPIRATION and BIOPSY paragraphs.

**Note:** In case of preceding aspiration, the biopsy specimen can be taken through the same incision in the skin; for taking the specimen, however, a different extraction point should be selected.

**NOTICE:** The biopsy needles can be locked with a luer cap.

**Intended use, aspiration:**

- Adaptation of a syringe by means of a luer lock fitting (4) on the biopsy needle.
- Take an aspiration specimen by exerting negative pressure (withdrawing the syringe plunger).
- Remove the biopsy needle including the syringe from the puncture channel.
- Remove the syringe with the aspiration specimen from the biopsy needle.

**Intended use, biopsy (taking a biopsy specimen under negative pressure):**

- To drive the biopsy specimen into the biopsy needle, gently push the biopsy needle forward by exerting rotating pressure. Using the ejector (7) and its calibration markings, check whether sufficient usable specimen material is in the biopsy needle (repeat Point 9 if required).
- Adaptation of a syringe by means of a luer lock fitting on the biopsy needle.
- Fix the tissue specimen in the biopsy needle by exerting negative pressure (withdrawing the syringe plunger).
- Remove the biopsy needle including the syringe from the puncture channel.
- Remove the syringe from the biopsy needle.
- Using the ejector, push the tissue specimen out of the biopsy needle.

**General remarks:**

- The product must be kept away from light and stored in a dry and cool place.
- This bone marrow biopsy system is a sterile single-use product intended for use on one patient. Resterilisation and reuse are not permitted.
- The reuse of single-use products involves a potential risk of infection for the patient or operator. Contaminated products can cause damage to health, disease or even death of the patient. Cleaning, disinfection or sterilisation can adversely affect essential material properties and product parameters and thus result in product failure.
- The attending doctor is responsible for the proper use and execution of the examination using this product (these operating instructions only serve as a general guideline).
- Disposal takes place in accordance with the hospital ordinance or the applicable laws. Due to the risk of infection and injury, care must be taken to avoid all contact with the tips or sharp edges of the instruments. This is especially important for contaminated instruments.
- Möller Medical GmbH expressly disclaims any liability for damage caused through improper use of the products or use by unauthorised personnel.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (only for European Union). Other countries should follow local regulations or procedures.

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