
BOW-1076
Bone Marrow Aspiration Biopsy

Medizinprodukte | Medical devices

Hersteller | Manufacturer

Vertriebspartner | Sales partner

Bestellnummer | Catalogue number

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Verwendbar bis [JJJJ-MM-TT] | Use by [YYYY-MM-DD]

Herstellungsdatum [JJJJ-MM-TT] | Date of manufacture [YYYY-MM-DD]

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
BOW-1076 Bone Marrow Aspiration Biopsy

en Instruction for Use for BOW-1076 bone marrow aspiration biopsy

Indication for Use: The Möller Medical Biopsy Needles and Systems are used to obtain core samples and aspiration from tissue by biopsy. The various wall thicknesses, diameters, tip geometries and cutting techniques of the different hollow needle types and biopsy systems affect the quality of the biopsy specimen.

Indications for the BOW-1076: The BOW-1076 bone marrow aspiration biopsy set is intended for aspiration in the area of the sternum or iliac crest. Application areas of the products are bone respectively bone marrow.

Contraindications: The use of the product, in particular on patients with blood coagulation disorders or patients treated with anticoagulants, is at the discretion of the attending doctor.

Patient Target Group: The device is applicable to all patients who, at the discretion of the physician, may be treated in compliance with the intended use, indications, contraindications, warnings and precautions.

Warnings and Precautions:

- The device is exclusively for use by authorised, specialised doctors who are familiar with the potential side effects, typical diagnostic findings, restrictions, indications and contra-indications of bone marrow aspiration biopsies.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (only for European Union). Other countries should follow local regulations or procedures.

Potential Complications: Potential complications associated with aspiration and core biopsy procedures are site specific and include, but are not limited to: cell seeding, infection, pain, bleeding, hematoma, non-target tissue or adjacent tissue injury, and organ or vessel perforation.

Preparatory Measures: Take the product out of the package and then remove the protective cover. Check the stylet tip and the biopsy cannula's cutting edge for damage and make sure that the stylet (1) slides freely within the biopsy cannula (2). **NOTE:** The arrow on the stylet cap and the cutting edge of the stylet tip, and also the arrow on the handle and the cutting edge of the biopsy cannula, point into the same direction. **CAUTION:** The product must not be used in case of damage or any other defect which can adversely affect its intended use.

Mode of Use:

1. Position the patient appropriately.
2. Definition of the entry point and marking if required.
3. To prevent infection, make sure to work under sterile conditions, including skin disinfection, sterile covering and sterile gloves.
4. Local anaesthesia of the puncture point.
5. If required, perform an incision using the scalpel in order to facilitate skin penetration.
6. Hold the biopsy needle (biopsy cannula with assembled stylet) by the handle, place your index finger on the needle and extend the finger towards the distal end. This enables safe handling and control over the maximum penetration depth. **NOTE:** Penetration depth is steplessly variable using the spacer ring (3), e.g. for sternal application; spacer ring may be removed, e.g. for iliac crest puncture.
7. Using the biopsy needle, advance forward through the incision and slowly push it with clockwise and counterclockwise rotary motion to form a puncture channel which first extends to the periosteum (resistance is noticeable), and then continue with the utmost care until you have overcome the compacta (the resistance suddenly decreases). **CAUTION:** Do not exert excessive force during application in order to prevent deformation or breakage of the biopsy needle causing inadvertent injury to the patient.
8. After entering the bone marrow cavity, remove the stylet from the biopsy cannula by rotating the locking mechanism 90°.
9. Adaptation of a standard luer lock syringe (not contained in the kit) by means of a luer lock fitting on the biopsy cannula.
10. Take an aspiration specimen by exerting negative pressure (withdrawing the syringe plunger).
11. Remove the syringe with the aspiration specimen from the biopsy cannula.
12. When the procedure is complete, remove the biopsy cannula from the puncture channel.

General Remarks:

- The device is supplied sterile. Before opening, check that the packaging is intact. Remove items from the packaging while maintaining aseptic conditions. Do not use if the packaging seal is not completely closed or if the sterile barrier has been damaged.
- The device should be stored in a clean, dry environment and protected from direct sunlight. Do not use after the expiration date on the label.
- The device is intended for single-use only and may not under any circumstances be reused. The reuse of a single-use device poses an infection risk to patients and users. Contamination of the device can lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization can adversely affect important material and design properties and device function.

- Disposal takes place in accordance with the hospital ordinance or the applicable laws. Due to the risk of infection and injury, care must be taken to avoid all contact with the tips or sharp edges of the instruments. This is especially important for contaminated instruments.
- The attending doctor is responsible for the proper use and execution of the examination using this product (these operating instructions only serve as a general guideline).
- Möller Medical GmbH expressly disclaims any liability for damage caused through improper use of the products or use by unauthorised personnel.