Möller Medical GmbH
Wasserkuppenstraße 29-31
36043 Fulda, Germany
Tel.: +49 (0) 661-9 41 95-0
Fax: +49 (0) 661-9 41 95-850
info@moeller-medical.com
www.moeller-medical.com

de

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en

Reprocessing

fr

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it

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es

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pt

Reprocessamento

nl

Herverwerking

da

Klargøring

sv

Rekonditionering

no

Reprosessering

fi

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tr

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el

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Möller Medical Blue

en Reprocessing of Möller Medical Blue

General fundamentals and references

Clean, disinfect and sterilize all devices before use. All devices are delivered non-sterile. (Cleaning and disinfection after removal of transport packaging, sterilization after packaging). Effective cleaning and disinfection has to be regarded as crucial for a valid sterilization.

Referring to your responsibility regarding sterility of devices in use ensure:

- Only validated and adequate procedures for cleaning/disinfection and sterilization of devices and products are to be used
- Employed devices (RDG, sterilizer) have to be maintained and checked regularly
- Validated parameters must be met each cycle

Strictly keep to instructed concentration, temperature, operating times and guidelines for rinsing regarding cleaning and disinfection media. Freshly made solutions, only sterile or low-germ water (max. 10 germ/ml) as well as endotoxin (max. 0,25 endotoxin-units/ml) water (e.g. purified water/ highly purified water) have to be used. For drying soft, clean and lint-free tissue (cloths) or filtered air has to be used.

Furthermore consider your country's legal regulations as well as hygienic rules for private clinics or hospitals. Especially this counts for different rules regarding an effective inactivation of prions (not applicable for the US).

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Cleaning and Disinfection

Basics

For an optimal cleaning and disinfection result an automatic process (RDG (cleaning and disinfection device)) has to be in charge. A manual procedure – even though ultrasound treatment is applied – should only be an exceptional case because effectiveness and reproducibility is clearly lower.

Pretreatment is mandatory in both cases.

Because of danger of splashing it has to be taken care that pretreatment is not made under running water.

Pretreatment

Remove rough contamination from devices immediately after indication (within max. 2 hours).

- Open and unlock devices (see specific instruction)
- Rinse devices for at least 1 minute in waterbath (temperature < 35 °C / 95 °F)
- Place the open and uncocked devices into preclearance bath 1 for the predefined operation-time and make sure that devices are sufficiently covered inside the bath. Check that devices are placed separately without touching each other.
- Assist preclearance by at least 3 times loading and triggering the devices while covered in the bath as well as by complete brushing of all inner and outer surfaces (each at the beginning of the predefined operation-time)
- Activate ultrasound afterwards (again for the minimum predefined operation-time, no less than 5 minutes)
- Then take devices out of the preclearance bath und rinse them at least 3 times properly (min. 1 minute) inside the waterbath

When choosing the detergent¹ please check that:

- It is generally suitable for the cleaning of devices made from metal and plastic
- The detergent is compatible with the ultrasound cleaning (no foaming)
- The detergent is compatible with the devices (see chapter "material resistance")

Machine Cleaning/Disinfection (RDG (Cleaning and Disinfection Device))

When choosing the RDG's Cleaning and Disinfection Devices please check that:

- RDG provides a proven effectiveness (for example DGHM- or FDA approval / Clearance / registration or CE sign acc. To DIN EN ISO 15883)
- Preferably a proven program for thermal disinfection (A0-value > 3000 or – if devices are older – min. 5 min at 90 °C / 194 °F) is in use (with chemical disinfection there is the risk of disinfection-agents residuals on the devices)
- The program is compatible with the devices and that it provides sufficient rinsing-cycles
- Only sterile or low-germ (max. 10 germs/ml) and low-endotoxin (max. 0,25 endotoxin units/ml) water (for example purified water / highly purified water) will be used
- Filtered air (oil-free, low-germ- and low-particle) will be used for drying and
- RDG will be maintained and checked regularly

When choosing the detergent-system please check that:

- It is generally suitable for the cleaning of devices made from metal and plastic
- There is an additional suitable disinfection-agent with proven effectiveness in use (for example VAH/DGHM- or FDA / EPA registration/clearance resp. CE sign) when there is no thermal disinfection applied. Suitable disinfection-agent has to be compatible with the detergents and
- That applied chemicals are compatible with the devices (see chapter "material resistance")

Procedure

- Open and uncock the devices (see specific instruction)
- Place the open and uncocked devices with the open side at the bottom inside the RDG. Check that devices are placed separately without touching each other.
- Start the program
- Remove devices from RDG after program end
- Check and pack devices as soon as possible after removal (see chapter "Check", "Maintenance", and "Packing") after an additional drying session in a clean area as appropriate

General proof of suitability for effective machine cleaning and disinfection of the devices has been provided by an independent, regulatory accredited and authorized (§15 (5) MPG) testing laboratory using RDG's G 7836 CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and pre-cleaning agent and detergent Neodisher MediZym (Dr. Weigert GmbH & Co. KG, Hamburg). The above described procedure was applied.

Manual Cleaning and Disinfection

When choosing the detergents and disinfection-agents please check that:

- They are generally suitable for the cleaning of devices made from metal and plastic
 - The detergent – if applicable – is compatible with the ultrasound cleaning (no foaming)
 - There is a suitable disinfection-agent with proven effectiveness in use (for example VAH / DGHM- or FDA / EPA registration/clearance resp. CE sign) and that this disinfection-agent is compatible with the detergents in use
 - That applied chemicals are compatible with the devices (see chapter "material resistance")
- Preferably combined detergents/disinfection-agents should not be used. Only in cases of minor contamination (no soiling visible) combined detergents/disinfection-agents can be used.

Cleaning

- Open and uncock the devices (see specific instruction)
- Place the open and uncocked devices into preclearance bath for the predefined operation-time and make sure that devices are sufficiently covered inside the bath. Check that devices are placed separately without touching each other. Assist preclearance by at least 3 times loading and triggering the devices while covered in the bath as well as by complete brushing of all inner and outer surfaces (each at the beginning of the predefined operation-time)
- Activate ultrasound afterwards (again for the minimum predefined operation-time, no less than 5 minutes)
- Then take devices out of the preclearance bath und rinse them at least 3 times properly (min. 1 minute) with water. Load and trigger the devices min 3 times inside the water.
- Check devices (see chapter "Check" and "Maintenance")

Disinfection

- Place the open, uncocked, cleansed and checked devices into disinfection bath for the predefined operation-time and make sure that devices are sufficiently covered inside the bath. Check that devices are placed separately without touching each other. Movable parts to be swayed several times during disinfection. Load and trigger the devices min 3 times inside the bath.
- Then take devices out of the disinfection bath und rinse them at least 5 times properly (min. 1 minute) with water. Load and trigger the devices min 3 times inside the bath.
- Dry devices by blowing of and out with compressed (filtered) air
- Pack devices as soon as possible after removal (see chapter "Packing") after an additional drying session in a clean area as appropriate

General proof of suitability for effective manual cleaning and disinfection of the devices has been provided by an independent, regulatory accredited and authorized (§15 (5) MPG) testing laboratory using pre-cleaning agent and detergent Cidezylme/Enzol and disinfection-agent Cidex OPA (Johnson & Johnson GmbH, Norderstedt). The above described procedure was applied.


Check

After Cleansing resp. cleansing and disinfection check all devices for corrosion, damaged surfaces, split-offs, contamination and discoloration. Select damaged devices (limitation to re-use see chapter "Re-use"). Devices still contaminated have to undergo cleansing and disinfection again.

Maintenance

When choosing device-oils please check that only those device-oils (white oil) that are approved for steam sterilization (considering the maximum applied sterilization temperature), are free from corrosion inhibitors and come with proven biocompatibility are used.

Parts to be oiled can be seen in picture 1:

Picture 1: parts to be oiled

Only apply the minimum quantity (spraying the complete device is forbidden!).
Device-greases are strictly forbidden!

Packing

Please pack devices resp. sterilization trays using only disposable sterilization packaging (single or double packaging) that comply with the following requirements (material/process):

- DIN EN ISO/AAMI ISO 11607 (for the US: FDA Clearance)
- Suitable for steam sterilization (temperature resistant to min. 142 °C (288 °F) sufficient steam permeability)
- Devices resp. sterilization packaging must have sufficient protection from mechanical damage

Sterilization

For sterilization only the following sterilization procedures are to be used; other sterilization procedures must not be used.

Steam Sterilization

- Fractional vacuum procedure^{2,3}, (with sufficient product drying⁴)
- Steam sterilization according to DIN EN 13060/DIN EN 285 resp. ANSI AAMI ST79 (for the US: FDA Clearance)
- Validation according to DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance evaluation (PQ))
- Maximum sterilization temperature 138 °C (285 °F, plus tolerance according to DIN EN ISO 17665)
- Sterilization time (exposure time with sterilization temperature)

country	fractional vacuum procedure	gravitation procedure
Europe, other countries	Min. 5 min ⁵ at 132 °C (270 °F) / 134 °C (273 °F)	non-recommended
USA	min. 4 min at 132 °C (270 °F), drying time min. 20 min ⁴	non-recommended

General proof of suitability for effective steam sterilization of the devices has been provided by an independent, regulatory accredited and authorized (§15 (5) MPG) testing laboratory using a steam sterilizer HST 6x6x6 (Zirbus technology GmbH, Bad Grund) and using a fractional vacuum procedure. Typical private clinic's and hospital's conditions were incorporated and the above described procedure was applied.

The "Flash sterilization" is forbidden

Please also do not apply hot-air sterilization, radio-sterilization, formaldehyde sterilization- or EO sterilization, as well as plasma sterilization

Storage

After sterilization all devices have to be stored in the sterilization package in a dry and dust-free area.

Material resistance

When choosing detergents and disinfection-agents please make sure to refrain from the following ingredients:

- Organic, minerally and oxidizing acids (minimum ph-value allowed 5,5)
- Alkaline solution (maximum ph-value allowed 8,5, neutral/encymatic cleaner recommended)
- Organic solvents (for example alcohol, ether, ketone, benzene)
- Oxidizing agents (for example hydrogen peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic/halogenated hydrocarbon

Never clean devices with a metallic brush or steel wool.

Rinse aid and acid neutralization agents must not be used.

All devices may only be exposed to temperatures below 142 °C (288 °F)

Re-use

Frequent reprocessing has little impact to these devices. The end of the product's operational life span is typically determined from attrition and use.

Fault removal

problem	Potential cause	Corrective actions
Biopsy device to be loaded with difficulty or not to be loaded at all	<ul style="list-style-type: none">• Incorrect handling• Clamping lever has not been pulled back completely• Safety button in wrong position• Device contaminated• Device is faulty	<ul style="list-style-type: none">• Correct handling• Clamping lever always to be fully pulled back• check safely button's position• clean and oil device• device to be sent back to manufacturer
Biopsy device cannot be triggered or can only be triggered with difficulties	<ul style="list-style-type: none">• device is not fully loaded• device is in locked position• biopsy device is contaminated• insufficient lubrication• device is faulty	<ul style="list-style-type: none">• check if device is fully loaded• check if device is in unlocked position• clean device• oil device• device to be sent back to manufacturer
Biopsy device does not fully release	<ul style="list-style-type: none">• Device is contaminated• Insufficient lubrication• Device faulty	<ul style="list-style-type: none">• Clean device• Oil device• device to be sent back to manufacturer

Möller Medical recommends to clean and lubricate the Möller Medical Blue after each usage!

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¹ If you – e.g. because of industrial safety reasons – use detergents and disinfectants media please note that they are free from aldehyde (otherwise fixation of blood contamination), come with proven efficiency (for example VAH/DGHM- or FDA/EPA registration/clearance resp. CE sign), suitable for disinfection and compatible with the devices (see chapter "material resistance"). Please note that the disinfectants in use during preclearance serve personal security; it does not replace disinfections to be done later on after cleansing took place.

² Min. three vacuum cycles

³ The use of the less effective gravitation procedure is only permitted when fractional vacuum procedure is not available and it requires generally far longer sterilization time as well as a product-, device-, procedure- and parameter-specific validation to be made under the solely responsibility of the user

⁴ The actual drying time needed is subject to the parameters that lie in the responsibility of the user (loading configuration and -density, sterilizer's condition, ...) and has therefore to be determined by the user. Nevertheless drying times should not fall below 20 minutes.

⁵ Resp. 18 minutes (prions'-inactivation, not applicable for the US)

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