# es

**REF**

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**QTY**

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Longitud total

Utilizar antes de AAAA-MM-DD

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**LOT**

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**STERILEEO**

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en

## Infusion Test Set *LiquoGuard*® 7

Safe for use in connection with  
CSF drainage pump *LiquoGuard*® 7

### ATTENTION:

Observe the *LiquoGuard*® 7 operating instructions.  
Do not use product if package is damaged.

**REF: 00 003 499**

Tube length: 1600 mm

### ATTENTION:

To use this tube set with the CSF drainage pump *LiquoGuard*® 7, please perform the following steps:

- Firstly, remove from its packaging a tubing set suitable for the lumbar infusion test and connect the sensor cable to the turned off *LiquoGuard*® 7.
- Restart the *LiquoGuard*® 7 via the On/Off switch.
- As soon as the tubing set is recognised by the *LiquoGuard*® 7, a green window will appear on the screen with the instruction to insert the tubing set.
- Select the button **Start Application** in the message window. The device is now in lumbar infusion test mode.

The *LiquoGuard*® 7 supports the operator in further preparatory procedures. Each measure is indicated in the Info bar of the display. Use the arrow keys in the dialogue window DIAGRAM to follow the installation procedure step-by-step.

Once preparations are completed, proceed with measurement.

- After pressing on the softkey Pause/Start to activate the pump, the pressure increases until a new stable pressure plateau is reached. This constant pressure that corresponds to the set speed of the pump is required to calculate the ROF.
- Once it is certain that the plateau value is stable and will no longer increase, enter the plateau pressure in the *LiquoGuard*® 7. To do this, press the softkey Plateau and change the value using the + and - buttons until it shows the same value as the stable Pcsf.
- The ROF for the application is indicated in the Info bar subject to Pstart and the plateau.

**Single-use:** Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning, disinfection and sterilisation may compromise essential material and design characteristics leading to device failure.

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