
# FAT Tubing Liposat<sup>®</sup> Pro plus

(REF. 00003948)

## Inhalt

de	Gebrauchsanweisung .....	3
en	Instructions for use.....	5
bg	Инструкция за употреба.....	7
cs	Brugsanvisning .....	9
da	Brugsvejledning .....	11
el	Οδηγίες χρήσης .....	13
es	Instrucciones de uso .....	15
et	Kasutusjuhend .....	17
fi	Käyttöohje .....	19
fr	Notice d'utilisation .....	21
hr	Upute za uporabu .....	23
hu	Használati utasítás.....	25
is	Notkunarleiðbeiningar .....	27
it	Istruzioni per l'uso .....	29
lt	Naudojimo instrukcija .....	31
lv	Lietošanas instrukcija.....	33
nl	Handleiding .....	35
no	Bruksanvisning.....	37
pl	Instrukcji obsługi .....	39
pt-BR	Instruções de uso.....	41
ro	Instrucţiuni de utilizare .....	43
sk	Návod na použitie .....	45
sl	Navodila za uporabo .....	47
sv	Bruksanvisning.....	49
tr	Kullanım Kılavuzu .....	51

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## en Instructions for use

The instruction for use is part of the product.

### Explanation of the safety symbols used

	Consult instructions for use		Do not re-sterilise
	Article number		Do not use if packaging is damaged
	Batch		Packaging unit
	Medical device		Use by YYYY-MM-DD
	Unique identifier of a medical device		Temperature limitation
	Sterilised using ethylene oxide		Store in a dry place
	Single sterile barrier system		Keep away from sunlight
	Single sterile barrier system with protective outer packaging		Date of manufacture YYYY-MM-DD
	Do not reuse		Manufacturer

**R<sub>X</sub> ONLY** Caution: Under US federal law, this device may only be sold to or ordered by a medical doctor.

Further information about the symbols used can be found on our website:

<https://moeller-medical.com/ifu-ga/ifu-ga/>

### Caution

- The product may only be used in combination with products from the Möller Medical aesthetics portfolio.
- These instructions are only valid together with the instructions for use of the respective devices used.
- Indications, contraindications, target group: See the instructions for use of the respective devices used.
- All incidents or serious adverse events occurring in relation to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is domiciled.
- Product cross-reactions are not known at this time.
- The attending doctor is responsible for the proper use and performance of the examination with this product.
- These instructions for use are only general instructions.
- Möller Medical GmbH accepts no liability for damage due to improper use of the products or the use of the products by unauthorised personnel.
- As the product is used in technically complex procedures, doctors should be familiar with the handling and use of the product and of the other instruments necessary for the procedure and the method.
- It is a sterile, single-use product.
- Check the integrity of the packaging.
- Contaminated products can damage health and cause illness or even the death of the patient.
- Do not use products with damaged, opened or moist packaging.
- The re-use of a single-use product presents a potential risk of infection for the patient or user.
- Cleaning, disinfection and sterilisation may compromise essential material properties and product parameters, leading to failure of the product.
- The products are to be disposed of in accordance with hospital directives/applicable laws.
- The product must be stored in a dry, cool, clean place protected from light.
- Check the expiry date.
- Do not use products after the expiry date.

## Intended use

The Liposat® Pro/power FAT tubing is used both for medical indications, including those accompanied by a change in fatty tissue, and for aesthetic body shaping. The Liposat® Pro plus FAT tubing is used to administer tumescent local anaesthesia, other aqueous infusion solutions, as well as endogenous subcutaneous tissue and its components into the body. The tubing set, sterilised with ethylene oxide, is marketed as a single-use product and may be used only in conjunction with the Liposat® Pro plus.

## Use and operation

1. Open the front cover of the pump.
2. The rotor rotates to a suitable position on its own.
3. Remove the sterile tube set from the packaging.
4. The red funnel of the tubing is on the left side.
5. Set the tubing length so that it can be connected to the infiltration container.
6. Press the tubing set into the left tube clamp until it clicks into place.
7. Insert the tubing so that it can be optimally picked up by the rotor.
8. Press and hold the "Insert tube set" button.
9. The rotor turns and takes up the tubing set.
10. Release the button when the tubing is threaded.  
Alternatives to 8-10:
  11. Turn the rotor to the right by hand.
  12. Carefully press the remaining tubing set into the right tube clamp until it clicks into place.
  13. Connect the red funnel to the outlet of the sterile infiltration container.
  14. Connect the right end of the tubing to the desired infiltration option.
15. Close the front cover of the pump.  
Left: Red funnel (infiltration medium)  
Right: Open end (infiltration option)

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