
TLA Tubing Liposat[®] Pro/power

(REF. 00002251)

Inhalt

de	Gebrauchsanweisung	3
en	Instructions for use.....	5
bg	Инструкция за употреба.....	6
cs	Brugsanvisning	9
da	Brugsvejledning	11
el	Οδηγίες χρήσης	13
es	Instrucciones de uso	14
et	Kasutusjuhend	17
fi	Käyttöohje	19
fr	Notice d'utilisation	20
hr	Upute za uporabu	23
hu	Használati utasítás.....	25
is	Notkunarleiðbeiningar	27
it	Istruzioni per l'uso	29
lt	Naudojimo instrukcija	31
lv	Lietošanas instrukcija.....	33
nl	Handleiding	35
no	Bruksanvisning.....	36
pl	Instrukcji obsługi	39
pt-BR	Instruções de uso.....	41
ro	Instrucţiuni de utilizare	43
sk	Návod na použitie	45
sl	Navodila za uporabo	47
sv	Bruksanvisning.....	49
tr	Kullanım Kılavuzu	51

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en Instructions for use

The instruction for use is part of the product.

Explanation of the safety symbols used

	Consult instructions for use		Do not re-sterilise
	Article number		Do not use if packaging is damaged
	Batch		Packaging unit
	Medical device		Use by YYYY-MM-DD
	Unique identifier of a medical device		Temperature limitation
	Sterilised using ethylene oxide		Store in a dry place
	Single sterile barrier system		Keep away from sunlight
	Single sterile barrier system with protective outer packaging		Date of manufacture YYYY-MM-DD
	Do not reuse		Manufacturer

R_X ONLY Caution: Under US federal law, this device may only be sold to or ordered by a medical doctor.

Further information about the symbols used can be found on our website:

<https://moeller-medical.com/ifu-ga/ifu-ga/>

Caution

- The product may only be used in combination with products from the Möller Medical aesthetics portfolio.
- These instructions are only valid together with the instructions for use of the respective devices used.
- Indications, contraindications, target group: See the instructions for use of the respective devices used.
- All incidents or serious adverse events occurring in relation to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is domiciled.
- Product cross-reactions are not known at this time.
- The attending doctor is responsible for the proper use and performance of the examination with this product.
- These instructions for use are only general instructions.
- Möller Medical GmbH accepts no liability for damage due to improper use of the products or the use of the products by unauthorised personnel.
- As the product is used in technically complex procedures, doctors should be familiar with the handling and use of the product and of the other instruments necessary for the procedure and the method.
- It is a sterile, single-use product.
- Check the integrity of the packaging.
- Contaminated products can damage health and cause illness or even the death of the patient.
- Do not use products with damaged, opened or moist packaging.
- The re-use of a single-use product presents a potential risk of infection for the patient or user.
- Cleaning, disinfection and sterilisation may compromise essential material properties and product parameters, leading to failure of the product.
- The products are to be disposed of in accordance with hospital directives/applicable laws.
- The product must be stored in a dry, cool, clean place protected from light.
- Check the expiry date.
- Do not use products after the expiry date.

Intended use

The Liposat® Pro/power TLA tubing is used both for medical indications, including those accompanied by a change in fatty tissue, and for aesthetic body shaping. The Liposat® Pro/power TLA tubing is used to administer tumescent local anaesthesia, other aqueous infusion solutions, as well as endogenous subcutaneous tissue and its components into the body. The tubing set, sterilised with ethylene oxide, is marketed as a single-use product and may only be used in conjunction with the Liposat® power and the Liposat® Pro.

Use and operation

Liposat® power

The pump runs anticlockwise.

1. Open the front cover of the pump.
2. Turn the pump head to a suitable position.
3. Remove the sterile tubing set with turbid pump segment from the packaging.
4. The spike is on the right side.
5. Place the adapter of the suction side of the hose set in the right tube clamp.
6. Turn the rotor to the left.
7. The pump hose guides help to thread the hose.
8. Place the adapter of the suction side of the hose set in the left tube clamp.
9. Connect the left end of the hose set to the second hose set and an infiltration option.
10. Insert the spike of the tube set into the bag containing the infiltration medium.
11. Close the front cover of the pump.
 - Luer on left (infiltration option)
 - Spike on right (infiltration medium)

Liposat® Pro

The pump runs clockwise.

1. Open the front cover of the pump.
2. The rotor rotates to a suitable position on its own.
3. Remove the sterile tube set from the packaging.
4. The spike is on the left side.
5. Carefully press the left end of the pump segment into the left tube clamp until it clicks into place.
6. Insert the tubing so that it can be optimally picked up by the rotor.
7. Press and hold the "Insert tube set" button.
8. The rotor turns and takes up the tubing set.
9. Release the button when the tubing is threaded.
 - Alternatives to 7-9:
10. Turn the rotor to the right by hand.
11. Carefully press the remaining pump segment into the right-hand tube clamp until it clicks into place.
12. Connect the right end of the hose set to the second hose set and an infiltration option.
13. Insert the spike of the tube set into the bag containing the infiltration medium.
14. Close the front cover of the pump.
 - Spike on left (infiltration medium)
 - Luer on right (infiltration option)

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