














Bending instructions

for Möller Medical single use cannulas

This instruction on cannula bending applies only for the Möller Medical GmbH **Vibrasat**[®] (Pro and power) cannulas. Please take note of the **Vibrasat**[®] cannulas instructions for use.

 <p style="text-align: center;">  Vibrasat[®] Single-Use Cannulas ASC-SXXX-XX-XXXX </p>	 	
 <p> Vibrasat[®] Reusable Cannulas ASC-RXXX-XX-XXXX </p>	 	 
  		
 <p style="text-align: right;">Max 20°</p>		